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## Creating a Longitudinal View Report

Cognitrax Report by CNS Vital Signs					Test Date: October 28, 2013 11:20:03				
Subject Reference/ID: Sample					Administrator: Technician				
Age: 64					Language: English (United States)				
Total Test Time: 18:38 (min:secs)					Cognitrax Online Version 0.0.1				
Patient Profile:	Percentile Range				> 74	25 - 74	9 - 24	2 - 8	< 2
	Standard Score Range				> 109	90 - 109	80 - 89	70 - 79	< 70
Domain Scores	Subject Score	Standard Score	Percentile	VI**	Above	Average	Low Average	Low	Very Low
Verbal Memory	47	90	25	Yes		x			
Psychomotor Speed	138	91	27	Yes		x			
Processing Speed	39	94	34	Yes		x			
Simple Visual Attention	51	99	47	Yes		x			
Motor Speed	48	89	23	Yes			x		

Domain Dashboard: Above average domain scores indicate a standard score (SS) greater than 109 or a Percentile Rank (PR) greater than 74, indicating a high functioning test subject. Average is a SS 90-109 or PR 25-74, indicating normal function. Low Average is a SS 80-89 or PR 9-24 indicating a moderate level of deficit or impairment. Below Average is a SS 70-79 or PR 2-8, indicating a moderate level of deficit or impairment. Very Low is a SS less than 70 or a PR less than 2, indicating a severe level of deficit and impairment. Reaction times are in milliseconds. An \* denotes that "lower is better", otherwise higher scores are better. Subject Scores are raw scores calculated generated from the test data values of the individual subtests.

VI\*\* - Validity Indicator: Denotes a guideline for representing the possibility of an invalid test or domain score. "No" means a clinician should evaluate whether or not the test subject understood the test, put forth their best effort, or has a clinical condition requiring further evaluation.

Finger Tapping Test (FTT)	Score	Standard	Percentile	
Right Taps Average	51	94	34	The FTT is a test of tapping speed on each hand. The F scores indicate motor speed. Lower scores indicate motor speed are faster with their hands.
Left Taps Average	45	86	18	
Verbal Memory Test (VBM)	Score	Standard	Percentile	
Correct Hits - Immediate	14	115	84	Verbal Memory Test measures the ability to recognize and remember verbal stimuli in the presence of distractors. The test measures the ability to recognize, remember, and respond to "Correct Hits" refer to correct responses. Lower scores indicate impairment.
Correct Passes - Immediate	12	74	4	
Correct Hits - Delay	10	103	47	
Correct Passes - Delay	11	74	4	
Continuous Performance Test (CPT)	Score	Standard	Percentile	
Correct Responses	40	103	58	The CPT measures sustained attention or vigilance and choice reaction time. Most normal subjects obtain near-perfect scores on this test. A long response time may suggest cognitive slowing and/or impairment. More than 2 errors (total) may be clinically significant. More than 4 errors (total) indicate attentional dysfunction.
Omission Errors*	0	103	58	
Commission Errors*	1	97	42	
Choice Reaction Time Correct*	514	72	3	
Symbol Digit Coding (SDC)	Score	Standard	Percentile	
Correct Responses	42	96	40	The SDC test measures speed of processing and draw upon several cognitive processes simultaneously, such as visual scanning, visual perception, visual memory, and motor functions. Errors may be due to impulsive responding, misperception, or confusion.
Errors*	3	81	10	

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## Longitudinal View Graph Result Report Domain

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